



Small carbon fiber calibration kit

The calibration kit includes a 300 mm carbon fiber wand and an L-frame distributed in a cardboard box.

The 300 mm wand is ideal for small to medium sized volumes. Just like the larger 600 mm wand, it has spherical markers for maximum calibration accuracy.

Being carbon fiber, it does not easily bend or expand as temperature changes.

FEATURES

- Non-deformable carbon fiber wand
- Spherical markers for better accuracy
- Adjustable level
- Floor-level origo
- Wand size 300 mm
- Lightweight



Aluminum L-frame

The L-frame is robust and resistant.




Water level

Two water levels make it easy to get an upright coordinate system.




Level adjustment

Get the L leveled quickly with the level adjustment screws



Floor-aligned origo

Properly align the coordinate system with your lab floor.



Spherical markers

Get maximum calibration accuracy with spherical markers



Carbon fiber wand

Strong, yet lightweight, the 300 mm wand is ideal for small to medium sized volumes.

WHAT'S IN THE BOX

- L-frame
- Small wand (300mm)
- Calibration certificate
- Product information



QUALISYS



[qualisys.com](https://www.qualisys.com)

sales@qualisys.com

Qualisys Europe

Gothenburg, Sweden

Qualisys Americas

Chicago, USA

Qualisys Asia Pacific

Shanghai, China