



Wand and L-frame for 3D calibration

This kit includes a carbon fiber wand, a folding L-frame and an aluminum carry case. The wand is 600 mm with spherical markers for maximum calibration accuracy.

Being carbon fiber, it does not easily bend or expand as temperature changes. A wand replacement set is also available.

FEATURES

- Non-deformable carbon fiber wand
- Spherical markers for better accuracy
- Folding L-frame
- Adjustable level
- Floor-level origo
- Force plate alignment plates
- Wand size 600 mm
- Lightweight



Folding L-frame

The L-frame can be folded to save space.



Water level

Two water levels make it easy to get an upright coordinate system.



Force plate alignment

Quickly align the L with a force plate using the alignment plates.



Floor-aligned origo

Properly align the coordinate system with your lab floor.




Spherical markers

Get maximum calibration accuracy with spherical markers.



Level adjustment

Get the L leveled quickly with the level adjustment screws.



WHAT'S IN THE BOX

- L-frame
- Wand grip
- Wand 600 mm
- Aluminum case
- Calibration certificate



QUALISYS



qualisys.com

sales@qualisys.com

Qualisys Europe

Gothenburg, Sweden

Qualisys Americas

Chicago, USA

Qualisys Asia Pacific

Shanghai, China