



Passive Markers

Qualisys super-spherical markers are available in 10 different variants. The markers are designed for maximum roundness and with their solid build, they are both durable and lightweight for efficient subject application.

The markers are available in packs of 10, 50 and 100 pieces, including plastic bases. Prices on the reverse side of this document do not include taxes and shipping costs, which vary by location. To get a price including these costs, please contact your local distributor or sales@qualisys.com.

Replacement plastic bases are also available separately.

FEATURES

- Flat base option
- Available in a soft neoprene or hard plastic units
- Screw-on option
- Available in a variety of sizes
- Underwater (UW) markers for in water captures
- Lightweight

SOFT MARKERS

Soft markers are rather durable units. They consist of two major parts: the marker itself with a built-in screw and a two-sided neoprene base. The marker is a foam ball covered by a special scotch-lite reflective film. Each package consists of 10 markers.



SUPER-SPHERICAL MARKERS

Our markers are high quality, lightweight markers designed for maximum roundness and an M3/M4 center-threaded cylinder resulting in greater spherical accuracy and identification of the marker centroid. The markers are sold in packages with 10, 50 and 100 markers and include plastic bases. UW are suitable for under water.



ATTACHING THE MARKERS

The passive reflective markers are usually applied directly on the skin with double adhesive tape. If the user would like to create his or her own markers, there is a retro-reflective marker material available in form of a tape.

SPECIFICATIONS



								Underwater markers		
Size	19 mm	16 mm	14 mm	12.5 mm	9.5 mm	8 mm	6.5 mm	19 mm	16 mm	14 mm
Weight	2.5g	2.0g	1.5g	1.0g	0.5g	0.3g	0.2g	2.5g	2.0g	1.5g
Thread	M4	M4	M4	M4	M3	M3	M3	M4	M4	M4

The markers are available in packages of 10, 50 or 100 pcs. Read about our other markers on <https://www.qualisys.com/hardware/>